

Archiving


This software copies the results of the Omron board results database (including image data) to an external archive where it is not subject to automatic deletion.

The results data can be viewed with the aid of any modern Web Browser, and from there, exported as Excel or CSV files.

The copy process works automatically for boards tested/verified after the Archiver has been installed. Older data can be imported from the Omron database by pressing a single button. If, for any reason a board cannot be archived, the responsible AOI or REP station will be alerted via a pop-up message window.

BitTrace14 Archiver

Long term data storage for your Omron/CKD AOI/AXI and Paste systems.



The Archiver allows you to perform long-term storage of your test and verification data. The system automatically imports data from all configured systems for storage that covers years, including image data from AOIs and AXIs as well as 3D paste data. Access to the stored data is provided via any standard web browser. Additional exports are supported including XML, CSV and Excel, either automatically on a board-by-board basis, or via a manual command. Search your data and see results for a given barcode across all supported systems. The software is expandable via plug-ins, for example to provide customized statistics.

